

FAR DIAGNOSTICS

CATALOG 2008

CLINICAL CHEMISTRY KIT – ROUTINE

| PRODUCT | METHOD | CODE | KIT SIZE |
|---------------------------------|----------------------------------|------|-----------|
| Acid Phosphatase | Kinetic - Colorimetric | 7006 | 22x3 ml |
| Albumin liquid | BCG - Endpoint | 7008 | 6x125 ml |
| Alkaline Phosphatase | Kinetic - S.C.E | 7010 | 24x3 ml |
| Alkaline Phosphatase liquid | Kinetic - S.C.E | 7012 | 4x50 ml |
| ALT / GPT | Kinetic UV - IFCC | 7014 | 24x3 ml |
| | | 7016 | 3x25 ml |
| ALT /GPT liquid | Kinetic UV - IFCC | 7018 | 5x50 ml |
| Amylase CNPG3 liquid | Kinetic | 7026 | 24x3 ml |
| | | 7028 | 6x50 ml |
| | | 7029 | 6x25 ml |
| Amylase EPS | Kinetic - EPS blocked | 7022 | 24x3 ml |
| | | 7024 | 5x30 ml |
| Amylase EPS liquid | Kinetic - EPS blocked | 7146 | 4x50 ml |
| AST /GOT | Kinetic UV - IFCC | 7032 | 24x3 ml |
| | | 7034 | 3x25 ml |
| AST /GOT liquid | Kinetic UV - IFCC | 7036 | 5x50 ml |
| Direct Bilirubin liquid | Jendrassik modified | 7038 | 5x125 ml |
| Total Bilirubin liquid | Jendrassik modified | 7040 | 5x125 ml |
| Direct + Total Bilirubin liquid | Jendrassik modified | 7042 | 4x125 ml |
| Calcium MTB liquid | MTB- Endpoint | 7044 | 6x50 ml |
| Chlorides liquid | Direct thiocyanate | 7048 | 6x50 ml |
| Cholesterol liquid | Trinder - Endpoint | 7050 | 6x50 ml |
| | | 7051 | 6x25 ml |
| | | 7052 | 6x125 ml |
| Cholesterol HDL liquid | Direct - Selective detergent | 7054 | 4x25 ml |
| Cholesterol HDL PEG liquid | Precipitation - PEG 6000 | 7056 | 6x25 ml |
| Cholesterol LDL liquid | Direct - Selective detergent | 7058 | 4x25 ml |
| Cholinesterase | Trinder - Kinetic benzoylcholine | 7060 | 24x3 ml |
| Cholinesterase liquid | Kinetic Colorimetric | 7062 | 5x25 ml |
| CK-MB | Kinetic UV - Immuno inhibition | 7064 | 24x2.5 ml |
| CK-NAC | Kinetic UV - NAC activated | 7066 | 24x2.5 ml |
| | | 7068 | 5x20 ml |
| Copper Serum | 3,5 - Di Br PAESA no deprotein. | 7070 | 5x25 ml |
| Copper Urine | 3,5 - Di Br PAESA | 7072 | 3x25 ml |
| Creatinine liquid | Jaffè - Fixed time | 7074 | 6x125 ml |
| | | 7075 | 6x50 ml |
| Fructosamine | Nitroblue Tetrezolium | 7080 | 5x20 ml |
| Gamma-GT | Colorimetric - Kinetic | 7076 | 24x3 ml |
| Gamma-GT Lio | | 7077 | 5x40 ml |

| PRODUCT | METHOD | CODE | KIT SIZE |
|---|---------------------------------|-------------|-----------------|
| Gamma-GT liquid | Carboxy Glupa | 7078 | 4x50 ml |
| Glucose liquid | Trinder - Endpoint | 7082 | 6x50 ml |
| | | 7084 | 6x125 ml |
| Alfa-HBDH | Kinetic UV | 7086 | 24x3 ml |
| Iron Ferene | Endpoint - no deproteinization | 7090 | 5x50 ml |
| Iron nitro Paps | Endpoint - no deproteinization | 7092 | 5x50 ml |
| LDH | Kinetic UV | 7094 | 24x3 ml |
| LDH liquid | Kinetic UV | 7096 | 5x50 ml |
| Lipase liquid | Colorimetric | 7098 | 4x25 ml |
| | | 7100 | 4x50 ml |
| Total Lipids liquid | Colorimetric | 7102 | 6x125 ml |
| Magnesium Calmagite liquid | Calmagite - Endpoint | 7108 | 4x50 ml |
| Magnesium Xylidyl Blue liquid | Xylidyl Blue - Endpoint | 7110 | 6x25 ml |
| Phospholipids | Trinder - Endpoint | 7116 | 3x25 ml |
| Inorganic Phosphorous | Molybdenum blue - Endpoint | 7088 | 4x25 ml |
| Total Protein serum liquid | Biuret - Endpoint | 7126 | 6x50 ml |
| | | 7128 | 6x125 ml |
| Total Protein Coomassie urine liq. | Coomassie blue | 7130 | 6x50 ml |
| Total Protein Pyrogallol red urine liq. | Pyrogallol red | 7132 | 6x50 ml |
| TIBC | Precipitation - Magnesium carb. | 7134 | 4x25 ml |
| Triglycerides liquid | Trinder - Endpoint | 7136 | 6x50 ml |
| | | 7137 | 6x25 ml |
| | | 7138 | 6x125 ml |
| UREA Color | Berthelot modified - Endpoint | 7140 | 4x125 ml |
| UREA UV | Fixed time - UV | 7142 | 5x50 ml |
| UREA UV liquid | Kinetic UV | 7144 | 5x40 ml |
| | | 7145 | 5x20 ml |
| Uric Acid liquid | Trinder - Endpoint | 7147 | 6x25 ml |
| | | 7148 | 6x50 ml |
| Zinc | 5 Br-PAPS no deproteinization | 7154 | 4x25 ml |