

Crony Instruments S.r.l.
AIRONE 200 RA



AIRONE 200RA - Random Access Chemistry Analyzer. Fully computer controlled with smart X-Y-Z robotics and sophisticated software with graphic presentation done in windows, offers great versatility, efficiency and user friendly operation. The analyzer offers many special features including Liquid Sensors for both sample and reagent; it does not require the expensive disposable accessories and offers trouble free operation of 80 to 160 tests/hour, depending on method. Furthermore, it includes : - a new sophisticated software, that enables operation in both Random Access as well as in Smart Batch mode. Ideal for STATS, since every 15 to 20 minutes a number of completed Final Patient Reports are printed out. - 27 On Line Reagents of which 9 reagents are cooled and the remaining 18 are at room temperature - Open system with an unlimited number of user programmable methods - an extra fast Peltier system to accelerate incubation in flow-cell - and many more small yet important features are included.