



## NEPHSTAR PLUS



- Can test 3 samples simultaneously
- No need to do sample dilution
- No need to do calibration curve
- Automatic Mixing
- Temperature controlled reaction chamber
- Automatic blanking
- Automatic calibration
- 1000 results storage

3 assays available now

1. HbA1c
2. CRP whole blood
3. Urine-Microalbumin

### ***HbA1C assay***

One step latex-enhanced direct method, no need for determination of total hemoglobin

Results according to NGSP or IFCC

Time saving

Assay range: 3-15%

### ***CRP Assay in whole blood***

Assay range: 1-350 mg/L

Only 10-20 ul venous or capillary blood required

Results available with 2 minutes