

[Reagent kit for Nephstar Plus](#)



Article name	Code	Pack	Measuring Range		Availability
Nephstar Plus bCRP Kit	BK051	100T	1.0-350	mg/L	Available
Nephstar Plus Urine Microalbumin Kit	DK061	100T	10-220	mg/L	Available
NephstarPlus HbA1c Kit	DK071	50T	3-15	%	Available

HbA1C assay

One step latex-enhanced direct method, no need for determination of total hemoglobin

Results according to NGSP or IFCC

Time saving

Assay range: 3-15%

CRP Assay in whole blood

Assay range: 1-350 mg/L

Only 10-20 ul venous or capillary blood required

Results available with 2 minutes