

Uritest-50



Uritest-50 Urine Analyzer

Technical Data:

- Parameters: pH, Nitrite, Specific Gravity, Ascorbic Acid, Blood, Glucose, Bilirubin, Urobilinogen, Ketone, Leucocytes, Protein, Calcium, Creatinine, Microalbumin
- Principle: dual-wavelength reflectance photometry
- Two test modes: Normal mode, 60 samples/hour, Continuous mode, 120 samples/hour
- Memory capacity: 1000 results
- Strip: Uritest 10A/11A urine reagent strips
- Dimension: 290mm x 200mm x 100mm

Features & Functions:

- Cold light source facility with long life and high stability
- 12V DC power supply by the adapter
- Built-in thermal printer & port for an external printer
- Two test modes and speeds for selection
- Data transfer capability via RS232 port
- Barcode reader optional

Microalbumin:

- Forecast the diabetic nephropathy and early detection
- Reach the sensitivity of 50mg/L with the 85% accuracy
- The best supplement for the test of protein