



FULLY AUTOCLAVABLE

Micro Volume Pipette

- Fully Autoclavable
- Light Weight & Ergonomic Design
- Click stop volume setting digital system
- Calibration conforms to DIN 12650 & EN ISO 8655 standards ensuring high accuracy & precision



Accumax

FULLY AUTOCLAVABLE

Micro Volume Pipette

Variable Volume

Cat No	Volumes	Increments
AV-100	0.5 - 10 µl	0.1 µl
AV-200	5 - 50 µl	0.5 µl
AV-300	50 - 200 µl	1 µl
AV-400	200 - 1000 µl	5 µl
AV-500	10 - 100 µl	0.5 µl
AV-600	100 - 1000 µl	5 µl
AV-700	2 - 20 µl	0.1 µl
AV-800	20 - 200 µl	1 µl
AV-900	500 - 5000 µl	50 µl
AV-1000	1 - 10 ml	0.1 ml

Fix Volume

Cat No	Volume
AF-0005	5 µl
AF-0010	10 µl
AF-0020	20 µl
AF-0025	25 µl
AF-0050	50 µl
AF-0100	100 µl
AF-0200	200 µl
AF-0250	250 µl
AF-0500	500 µl
AF-1000	1000 µl
AF-2000	2000 µl
AF-5000	5000 µl
AF-10000	10000 µl



OEM
Option Available


ISO 9001 : 2000




C0508 ISO/IEC 17025:2005
National Accreditation Board for
Testing and Calibration Laboratories

ISO 17025:2005



www.accumaximum.com

Fine Care Biosystems